

DKW-320 Microplate



Main Features :

- * Large LCD display, menu designed for convenient use
- * Whole plate washing or single strip washing
- * 12-way and 8-way manifold included
- * Low residual volume by two pipettes
- * Complete bottom washing
- * Shaking and dipping function
- * Automatic monitoring of vacuum and pressure ,automatic rinse cycle
- * Flat, V-bottom or U-bottom plate and strips washing
- * Large memory to store up to 48 user programmed wash protocols
- * Emergency stop available and liquid warning

Technical Specifications :

Manifold: 8 pins and 12 pins
Wash time: 0~99 times adjustable
Wash strips: 1~12 adjustable
Dipping time: 0~3600s adjustable
Shaking time: 0~600s adjustable
Resolution: 50~3000ul/hole
Dispense precision: $CV \leq 2\%$
Sipping time: 0.1~9.9s
Power supply: 198~242V, 49~51Hz